



MONARK 827E



MAIN FEATURES

- Manual power control
- No power adaptor is needed

THE BASICS

- Good cycling feeling
- 5 kHz HR system (only in display)
- High quality components for heavy duty users
- Suitable for training and rehabilitation

Max user weight 250 kg

Medical Device class IIa

Specifications

MONARK 827E

Certifications	EU directive 93/42/EU class IIa
Frame	
Type	Monark Classic
Handlebar	Monark Handlebar
Stem	Monark Standard
Seatpost	Monark Standard
Saddle	Saddle, blue
BB	68/122 mm
Crank	Steel, 170 mm
Pedals	9/16", pedals with foot straps
RESISTANCE	
System	Monark Basic
Flywheel	20 kg
DISPLAY	
RPM	Yes
HR	5 kHz
Time	Yes
Speed	Yes
Distance	Yes
Watt	n/a
Connectivity	
RS232	n/a
Usb	n/a
WEIGHT & MEASUREMENTS	
Max user weight	250 kg
Length [mm]	1120 mm
Width [mm]	519 mm at handlebar 640 mm at support legs
Height [mm]	945-1295 mm at handlebar 780-1105 mm at saddle
Weight [kg]	47 kg
POWER REQUIREMENTS	Batteries in the display



HEART RATE

Equipped to measure heartrate via chestbelt.

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 15 36 01 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden